Glamatic Bipap - Breath Easy Auto-Bipap ~ Bi-Level Automatic Positive Airway Pressure

Model - GPBI107





Sleep Report on Screen



SD Card Slot for Data Storage & Report

Features:

2.8" Color TFT Display ensures a friendly user interface

One Touch Button operation makes the parameter setting easy.

Constant temperatures humidifire ensure comfortable airflow to patient.

ERP function automatically traces your breath, provides proper treatment to patient-

Continuous Sleep Rreport on Screen.

Latest Technology of Different pressure sensor & Automatic Leakages sensor makes the therapy more accurate

Customized German Motor and special sound Insulation material makes very low noise

Patient breathing with mask can trigger the machine easily into Working Mode:

Taking off the machine in working mode make it enter into Standby mode automatically 15 sec.

Internal Intelligence Flash records every using case including Real-time data, event information, cct.

Working Modes:

Auto BiPAP : CPAP, Auto CPAP, S, T & S/T

Pressure Range (IPAP) : CPAP - 4 to 20 cmH2O, BiPAP : 4 to 25 cmH2O Auto Pressure : Yes Automatic On/Off: Yes Ramp Time : 0-60 Min Auto Leackage Detection : Yes

Sound Level : <25dBa New Weight : 1.6 Kgs

: 25 x 14 x 9.5cm (Host + Humidifier) **Dimensions**

Package Includes:

- · Main Unit
- · Heated Humidifier
- · Nasal Mask
- · Tubing 6 Ft
- · Carry Case
- · Filter
- · SD Card for Unlimited Online Report











