Jlamatic

Oxygen Concentrator Model - GPOC109

Active lifestyle of your patients in today's world, demands easy and efficient solutions to their problems. Glamatic GPOC101 fits right in. It combines comfort and affordability to prove convenient oxygen therapy to your patients at home.

Oxygen therapy becomes easy, reliable and economic with Glamatic Concentrator over the course of its usage. To ensure comfortable experience during its ownership, it is incorporated with essentials like strong plastic case, proven sieve bed, low power consumption capacity and Pressure Swing Adsorption (PSA) technology. The compact design, light weight and advanced features make it ideal for convenient and hassle-free oxygen therapy for your patients at home.



High Purity of Medical Grade Oxygen High Grade Molecular Sieves to provide 93+% of Oxygen per liter



LCD Display

LCD Display shows the time the machine is run & various Alarms



Machine is equipped Automizer Therapy



NASA derived technology

Rapid Pressure Swing Absorption Technology



Adjustable Flowmeter

Flowmeter can be easily adjustable as per the need of patient from 0~5 Lpm



Triple Protection Filter

Equipped with Gross Particle Cabinet Filter, Compressor Intake Filter & Bacterial Filter



User Friendly / Zero Maintainence

Rugged & Durable ABS moulded cabinet for strength and durability



Low Cost of Oxygen Generation

No other raw material required Except Air



Audio & Visual Alarm

Equipped with Audio & Visual Alarm System



Low Noise

Quiet Unit-allowing restful sleep to patients.



Real-time Flowmeter Knob

Flow Control Kob allows you the ability to control prescribed oxygen level. High Output Allows longer tubing without restricted the flow.



Low Power Consumption

Consumes very less power











