

### SINGLE SURFACE WITH TROLLEY UST



#### Introduction :

18 Pc. Hi power LED provides to cover the bed area and for effective Phototherapy treatment. No ultraviolet rays ensures minimize the UV exposures. No infrared rays ensures minimize the water loss. Adjustable Height - Conveniently allows you to control Radiance levels on the neonate and different types of cradle. Convenient Tilting of Source Unit, Height Adjustment with Single hand and Castor with elegant maneuvering helps to place the system at different positions, at user convenience. Lamp usages hours greater than 45,000 hours ensures no need of frequent change of lamp.

#### SPECIFICATIONS

- ▶ Time Totalizer (Optional).
- ▶ Power LED's emitting blue light in 450 - 470 NM.
- ▶ Long Life > 45,000 hrs.
- ▶ Intensity: >25  $\mu\text{w} / \text{cm}^2 / \text{HM}$
- ▶ M.S. Fabrication with powder coated paint.
- ▶ Operation Voltage: 100-250 V AC, 50Hz.

### DOUBLE SURFACE WITH TROLLEY DSP



#### Introduction :

18 Pc. Hi power LED provides to cover the bed area and for effective Phototherapy treatment. No ultraviolet rays ensures minimize the UV exposures. No infrared rays ensures minimize the water loss. Adjustable Height - Conveniently allows you to control Radiance levels on the neonate and different types of cradle. Convenient Tilting of Source Unit, Height Adjustment with Single hand and Castor with elegant maneuvering helps to place the system at different positions, at user convenience. Lamp usages hours greater than 45,000 hours ensures no need of frequent change of lamp.

#### SPECIFICATIONS

- ▶ Time Totalizer (Optional).
- ▶ Power LED's emitting blue light in 450 - 470 NM.
- ▶ Long Life > 45,000 hrs.
- ▶ Intensity: >25  $\mu\text{w} / \text{cm}^2 / \text{HM}$ .
- ▶ M.S. Fabrication with powder coated paint.
- ▶ Operation Voltage: 100-250 V AC, 50Hz.